Patella Tracking Orthosis (PTO) Knee Brace

PTO-KNIESCHIENE
SUPPORTO PER GINOCCHIO PER L'ORTOSI PATELLARE (PTO)
ORTHÈSE DE GENOU POUR ALIGNEMENT ROTULIEN (PTO)
RODILLERA PARA ORTOSIS DE RÓTULA (PTO)




EC REP

E/U authorized representative
MDSS GmbH
Schiffgraben 41
D-30175 Hannover
Germany

Breg, Inc.
2885 Loker Ave. East
Carlsbad, CA 92010 U.S.A.
P: 800-321-0607
F: 800-329-2734
www.breg.com

AW-1.08770 REV G 0213

W A R N I N G

CAUTION: FOR SINGLE PATIENT USE ONLY.

CAUTION: IF YOU EXPERIENCE INCREASED PAIN, SWELLING, OR ANY ADVERSE REACTIONS WHILE USING THIS PRODUCT, IMMEDIATELY CONSULT YOUR MEDICAL PROFESSIONAL.

WARNING: CAREFULLY READ USE INSTRUCTIONS AND WARNINGS PRIOR TO USE.

WARNING: THIS DEVICE WILL NOT PREVENT OR REDUCE ALL INJURIES. PROPER REHABILITATION AND ACTIVITY MODIFICATION ARE ALSO AN ESSENTIAL PART OF A SAFE TREATMENT PROGRAM.

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A B R A C E A P P L I C A T I O N

Unfasten and loosen all straps and fasten hook/loop tabs back against the strap as shown. This will facilitate pulling on the brace.

While seated, pull on brace. The front hole must be positioned directly over the kneecap (patella), as shown. The plastic pressure plate should be positioned on the outside of the knee.

T I G H T E N F O U R M A I N S T R A P S

1st: lower rear; 2nd: upper rear;
3rd: lower front; 4th: upper front.

T I G H T E N I N G P A T E L L A B U T T R E S S S T R A P S

Bend the knee to 60°. While pushing down on the pressure plate, tighten the buttress straps around the kneecap. These straps should be snug.

NOTE: over-tightening will cause the medial (inside) hinge to pull forward.

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