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
1 • DESCRIPTION

The KINETEC 6080 is used for passive mobilization of the Elbow, providing the following motion:

- Flexion-extension of the Elbow from 0° to 135°.
- Simultaneous 90° pronation and 90° supination with flexion-extension.

A-Upper arm support • B-Forearm support • C-Telescopic column • D-Motor • E-Swiveling handle • F-Control box • G-Base.

2 • ASSEMBLY

- Attach the telescopic column (C) to the red base (G) with the lock nut and the washer.
- Place the wire basket (1) over the telescopic column (C).
- Attach motor (D) to the telescopic column with the toothed link (2).
- Place the control box (F) on the wire basket (1).

3 • ELECTRICAL CONNECTION: SAFETY FIRST.

The KINETEC 6080 is a Type B Class I machine.

Before plugging in the machine, make certain that the available voltage corresponds to the voltage indicated on the back of the control box (110 or 240 Vac 50/60Hz).

Connect the cords to the back of the control box as follows:

- Attach the hand control cord (3) to terminal (4).
- Attach the cord (5) from the the control box (6) to the motor (D).
- Attach the power cord (7) to the control box (8).

IMPORTANT:

Cords (3) and (5) must be securely fastened.

Check to be sure that the plug fits the socket being used.

The socket must be grounded.

Use only the power cord delivered with the KINETEC machine.

Make sure that the power cords are not tangled during operation.

Fuse change:

Warning: Replace the fuse (9) with one of the same type and value (see value indicated on the control box (F)).

4 • SAFETY

The healthcare professional determines the CPM protocol to be used (set up, duration and frequency of each session).

The patient must understand the Start/Stop/Reverse function of the hand control.

The KINETEC 6080 unit is CE marked.

DANGER: Risk of explosion if used with flammable anesthetics.

5 • START UP

Turn on the power switch (10). The green light will turn on.

START / STOP / REVERSE function

The KINETEC 6080 is equipped with a START/STOP/REVERSE function:

- Press the button once, the machine starts.
- Press the button twice, the machine stops
- Press the button a third time, the machine reverses direction..

IMPORTANT: for maximum safety, the hand control should always be given to the patient.

Adjusting the SPEED

The SPEED is set by turning the button (12). Minimum speed is all the way to the left; maximum speed is all the way to the right. The speed can be modified while the machine is running or stationary.

Adjusting the LOAD

The LOAD allows the physician to adjust the motor power to patient's need.

The LOAD is set by turning button (27).

Modification of the load can be done while the machine is running or stationary.

6 • ADJUSTING THE LIMITS

The physician positions the cursors (13) and (14) on the threaded rod and locks them with the safety nuts (15) and (16).

The graduated scale placed on the screw drive casing gives only an approximate measurement.

7 • USE OF THE KINETEC HYGIENIC PADS

Place the pads (17) and (18) on the adhesive parts of the forearm and arm support, with the sponge side up.

FOR OPTIMAL HYGIENE, A NEW SET OF PADS SHOULD BE USED FOR EACH PATIENT.

The straps (19) and (20) are used for holding patient's arm on the machine. Do not use strap (19) when machine is used for pronation-supination.

For glove use, see the drawing.

For cleaning, see paragraph 12.

Part number to order pads : 4650000646

Part number to order glove : 4650000654

8 • SETTING UP THE PATIENT

Make sure that clean pads are on the machine (see section 7 and 12).

Place the patient comfortably in a sitting or reclining position.

Place arm into the machine, in the maximum tolerable flexion.

Adjust the height of the machine on the telescopic column (C) and the incline of the toothed link (2).

Adjust the joint axis of elbow and the joint axle of the machine.

Lock casters of base (G).

Place the pronation-supination hand grip (E) in the palm of hand with hand in supine position.

Fix the top of arm with the strap (20)(Fig.7)

9 • POSSIBLE MOVEMENTS

Flexion-Extension of the elbow (H).

The KINETEC 6080 allows flexion-extension of the elbow from 0° to 135° with fixed PRONO-SUPINATION (H).

Choose the hand position between full PRONATION to 90° and full SUPINATION to 90°.

- Position the cleat (25) on the LEFT of its slide (Fig.9•1).
- Lock the cleat (25) with the locking nuts (21 and 22; Fig.9•1).
- Wind the rope according to the indicated direction on the primary pulley (23), and put the rope into the flexion pulley (24) and into the cleat (Fig. 9•2 or 9•3).
- Tighten the rope to bring the hand to the chosen position.
- Cinch the rope into the cleat (25).

Fig.9•2: representation for a left arm; Fig.9•3: representation for a right arm.

AT THE TIME OF FLEXION-EXTENSION, THE HAND-GRIP WILL STAY AT THE INITIAL ANGLE.

PRONATION-SUPINATION in conjunction with Flexion-Extension of the elbow (H+I).

The KINETEC 6080 allows PRONO-SUPINATION 90° to 90° in conjunction with flexion-extension of the elbow from 0° to 135° (H+I).

- Position the cleat (25) on the RIGHT of its slide (Fig.9•4).
- To obtain pronation-supination associated less than 180°, move the cleat to the left on the slide.
- The more the cleat is moved to the left, the less pronation-supination occurs.

Fig.9•5: representation for a left arm; Fig.9•6: representation for a right arm.

10 • MAINTENANCE

Direct a jet of compressed air into the protection screw casing.

Do not clean the screw with cotton. Move plug that holds up the nut, and grease the screw.

The motor does not need maintenance. However, compressed air can be used to remove carbon buildups. After 4,000 hours of operation, returning the machine to Smith & Nephew for maintenance is recommended.

11 • BREAKDOWN

If you need a spare parts catalog and/or service manual, please request it from your KINETEC distributor.

- After connecting the power cord on the KINETEC 6080 machine and switching ON the machine (10), if:

- The green light is off,
 - Try another electrical device to make sure that the power supply is uninterrupted.
 - Check the fuse of the socket (9) by replacing it with a fuse of the same type and value.
 - If the green light is still off, contact your KINETEC specialist.
- The KINETEC 6080 fails to operate but the green light remains on,
 - Press the START/STOP switch two more times.
 - If your KINETEC 6080 still fails to operate, contact your KINETEC specialist.
- The red light (26) is on,
 - Verify that the hand control is firmly connected. If so, contact your KINETEC specialist.

12 • CLEANING

First, switch the CPM unit OFF and remove the power cord from the receptacle.

Spray a DISINFECTANT product (e.g., PROPANOL/ISOPROPANOL or ALDEHYDE solution) on to the SURFACES (plastic support, plastic covers and metal parts).

The KINETEC HYGIENIC pads are machine-washable (30°C, 90°F) with chlorine bleach.

For optimal hygiene, a new set of pads should be used for each patient.

13 • SPECIFICATIONS

Product

Weight: 21 Kg (46.2 lbs)

Dimensions: 72cm x 65cm x 130cm
(25.3 x 25.6 x 51.2 inches)

Range of motion: from 0° to 135° of flexion with
90° pronation to 90° supination.

Speed: from 3'50" to 1'00" for a cycle from 0° to 135°

Patient sizing: from 1,40m to 1,90m
(4.5 to 6.25 feet)

Electrical

Power supply: 110/230 V~

Frequency : 50/60Hz

Electrical power : 50VA

Current consumption :
220 mA (under 230 V~)

110 mA (under 110 V~)

Type B Class I device

Environmental

- Storage / shipping conditions:

Ambient temperature: -40°C to +70°C
(-40°F to +160°F).

Relative humidity up to 90%.

- Operating conditions:

Ambient temperature: +10°C to +40°C
(+50°F to +105°F).

Relative humidity 80%.

14 • SYMBOLS

	TYPE B EQUIPMENT (insulated against electric shock)
	WARNING (consult accompanying documents)
	OFF (disconnected from the mains)
	ON (connected to the mains)

	Remote control switch
LOAD	LOAD
	Minimum SPEED
	Maximum SPEED

15 • GUARANTEE LIMITS

The KINETEC guarantee is strictly limited to free replacement or to repair of the part or parts acknowledged as defective, at the Smith & Nephew factory. Smith & Nephew Kinetec guarantees its passive motion device against any manufacturing defect for 1 year from the date of purchase by the consumer.

Smith & Nephew Kinetec alone is qualified to judge whether the guarantee applies to the specific problem.

The guarantee cannot be applied if the apparatus has been subjected to abnormal use or has been used in circumstances other than those stipulated in our user manual. The guarantee will not be applied in case of deterioration or accident due to negligence, lack of supervision or maintenance or an attempt to repair the equipment.

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