

# SUMMARY

Page

<b>DESCRIPTION</b>	<b>2-6</b>
<b>ASSEMBLY &amp; TRANSPORT</b>	<b>7</b>
<b>USE OF THE HAND CONTROL</b>	<b>8-18</b>
<b>INSTRUCTIONS FOR USE</b>	<b>19-31</b>
<b>ACCESSORIES</b>	<b>32</b>
<b>PRODUCT INFORMATION</b>	<b>33-35</b>

GB

## SAFETY RECOMMENDATIONS

Select your language  
SEE PAGE 36

- Before use, please read this document.
- The physician defines the protocol and ensures that it is correctly implemented (adjustments, session durations and frequency of use).
- We recommend that you lock the hand control when you give it to the patient.
- For optimum safety, always give the hand control to the patient before starting the system.
- Explosion hazard: KINETEC Centura is not designed for use in the presence of flammable anesthetics.
- Check that the electrical socket is in good condition and is suitable for the splint power supply cord. The latter complies with current standards and has a grounding socket. The plug may be connected to any standard socket. The socket must however have a grounded pin.  
To connect the power supply, only use the original cable supplied with the machine.  
Check that the cables remain free around the device so that they do not get damaged.

Manual N°: 467896270 – Updated April 2002

KINETEC Centura – Series 2

KINETEC and Centura are trademark of AbilityOne.



## DESCRIPTION

### DEFINITION

The KINETEC Centura is a upper extremity PASSIVE mobilization device enabling the following movements:

- Extension 20°                      Flexion 180°.
- Adduction 20°                      Abduction 160°.
- Internal Rotation 60°              External Rotation 90°.
- Synchronized Abduction + Rotation
  - Adduction 20°              Abduction 160°.
  - Internal Rotation 30°              External Rotation 90°.
- Horizontal Abduction from -30° to 110° (available on Centura 5)

### • Indications

- Total shoulder replacement.
- Repeated dislocation of the humerus.
- Rotator cuff injury.
- Upper humerus fractures
- Scapula fractures.
- Acromioplasty.
- Capsulotomy.
- Arthrolysis.
- Synovectomy for Rheumatoid Arthritis.
- All type of shoulder stiffness joint.

### • Clinical Benefits

- Breaks the cycle of trauma, inflammation and the loss of range of motion.
- Prevents joint stiffness.
- Speeds the recovery of post-operative range of motion.
- Maintains the quality of the articular surface.
- Reduces pain and edema.
- Promotes articular cartilage healing.
- Reduces hospitalization time
- Reduces the need for pain medication.

### • Contraindications


- Unstable fractures.
- Spastic paralyses.
- Uncontrolled infection.
- The machine are not adapted for patients height more 2m(6'7") or under 1.4m(4'7")

## DESCRIPTION

### ELECTRICAL CONNECTION

KINETEC Centura is a type B class I device.  
Before connecting the device to the power supply, check that the mains voltage matches that shown on the plate (100-240 V~ / 50-60 Hz) below switch ON (2).

Connect the power supply cable (1).



GB

### IMPORTANT



Check that the electrical socket is in good condition and is suitable for the splint power supply cord. The latter complies with current standards and has a grounding socket. The plug may be connected to any standard socket. The socket must however have a grounded pin.  
To connect the power supply, only use the original cable supplied with the machine.  
Check that the cables remain free around the device so that they do not get damaged.

The cables (motors and hand control) can be plugged in any of the connectors

### Starting the unit

Switch on (2).

While the unit begins an auto diagnostic, the display shows the following:



Your KINETEC Centura is ready to be used.

### SAFETY

The physician defines the protocol and ensures that it is correctly implemented (adjustments, session durations and frequency of use).

The patient must know the start/stop/reverse function on the control handle. Hand control must be accessible to patient at all times. See page 8.

KINETEC Centura complies with Directive 93/42/CEE and UL 2601.

### EXPLOSION HAZARD:

**KINETEC Centura is not designed for use in the presence of flammable anesthetics.**

In case of electromagnetic interference with other devices move the device.  
KINETEC Centura is in compliance with standards in force (IEC 601.1.2),  
electromagnetic compatibility standard for medical devices.



## DESCRIPTION

### DESCRIPTION • COMPONENTS

GB

KINETEC Centura consists of the following components:

- 1 – Chair
- 2 – Frame
- 3 – Wheels
- 4 – Transport handle
- 5 – Arm rest knobs
- 6 – Arm rest
- 7 – Hand control support




- 8 – Locking of the right/left sliding
- 9 – Locking of the up/down sliding
- 10 – Chair locking knob
- 11 – Locking of the abduction motor
- 12 – Abduction motor
- 13 – Shoulder depth sliding lock

**KINETEC**

AbilityOne


## DESCRIPTION

### Abduction associated with rotation splint




- 14 – Arm length setting lock
- 15 – 90° elbow splint
- 16 – 90° elbow splint lock
- 17 – Rotation motor
- 18 – Rotation motor lock
- 19 – Forearm length setting lock
- 20 – Forearm slider
- 21 – Right/left bean swivel lock
- 22 – Forearm splint

### Abduction or flexion splint



- 18 – Swiveling splint support lock
- 23 – Swiveling splint support
- 24 – Elbow flexion setting lock
- 25 – Arm splint
- 26 – Arm splint lock
- 27 – Forearm length setting lock

### Horizontal Abduction splint (available on Centura 5)



- 28 – Horizontal abduction column
- 29 – Arm length setting lock
- 30 – Elbow splint support lock
- 31 – Elbow support older
- 32 – 90° elbow splint lock
- 33 – Forearm length setting lock
- 34 – Forearm slider



GB

## DESCRIPTION

### DESCRIPTION • ELECTRICAL

GB

- 1 – Supply cable connector switch
- 2 – Fuse
- 3 – ON/OFF switch
- 4 – Hand control lock switch
- 5 – Defect or power light
- 6 – Motor or hand control connectors
- 7 – Hand control location for transport
- 8 – Hand control



- 9 – Liquid-crystal display
- 10 – Flexion/extension
- 11 – Abduction/adduction
- 12 – Rotation
- 13 – Abduction/adduction synchronized with rotation
- 14 – Lower limits setting
- 15 – Upper limits setting
- 16 – Increase / decrease
- 17 – START
- 18 – STOP
- 19 – FORCE
- 20 – SPEED
- 21 – PAUSE
- 22 – TIMER
- 23 – PROGRAM access
- 24 – Horizontal abduction

**KINEIEC**

AbilityOne

## ASSEMBLY & TRANSPORT




### Base assembly

Remove all the components from their packaging.

Unscrew the chair locking knob (3).

Position the chair (1) on the base (2) with the back of the chair at the wheel and screw the chair locking knob (3) to secure it in place.

Place the elevation motor (4) on right or left, depending on the limb be moved.



The other components to be used depend on the selected movement.



### Unit transport


For easy transport of the unit, it features 2 wheels (5) and a handle (6).



Place the arm support as close as possible to the chair to limit the overall dimensions and help balance the unit.

Place your foot as indicated to balance the unit.

You can adjust the height of the handle with knob (7).



## USE OF THE HAND CONTROL

### Locking the hand control setting

The hand control allows the patient to control the machine as appropriate.

The switch (4) has 3 positions:



#### LOCKED POSITION (1)

The operational settings can be read and the START/STOP/REVERSE function operated.





#### UNLOCKED POSITION (2)

All adjustments are possible.





#### HALF-LOCKED POSITION (4)

It is possible to switch the program and modify the upper and lower movement limits. The START/STOP/REVERSE function is always accessible.



#### Double locking

Simultaneously press the  and  keys to lock the hand control.



The display reads LOCK. To unlock the hand control, simultaneously press the same keys. The display reads UNLOCK.

You can not change the parameters, if you try the display reads:

LOCK 1: if locked with only the switch (4)

LOCK 2: if only double locked

LOCK 12: if locked with the switch (4) and double locked.

**We recommend that you lock the hand control when you give it to the patient.**

### START/STOP/REVERSE function

As with all KINETEC systems, KINETEC Centura is equipped with a START/STOP/REVERSE function.

When the unit is running, the display reads RUN

Press the



key of the hand control. The movement stops. The display reads STOP

Press the



key of the hand control. The movement starts in the opposite direction and the display reads RUN.

#### **Caution:**







**For optimum safety, always give the hand control to the patient before starting the system.**

## USE OF THE HAND CONTROL

### Reset time function



This function allows one to read the running time since the last resetting of the counter.

GB

Beginning	Keys to press	Display	Remarks
To stop the unit		FLEX/EXT 50 STOP 15 100	Check if the locking switch is in the following position  or 
Press simultaneously on the 2 keys	 	RESET TIME 125H Reset: limit low	The display indicates the running time since the last resetting.
To reset the counter, press the key		RESET TIME 125H Reset is done	The counter is now reset.
Or After 5 seconds, the reset function switches off and the running time remains in the memory.		FLEX/EXT 50 STOP 15 100	

### How to choose a movement








First switch the machine off  and put the switch in the  position in order to change the movement.

You can select a movement by pressing the appropriate button. The LED is on.  
**When a movement is first selected, the system returns to the original parameters of the movement (default settings).**

### Default settings:

Default settings:

						
	Abduction	Flexion/Extension	Rotation	Abduction + Rotation		Horizontal* abduction
• Lower limit	30°	30°	0°	30°	0°	-30°
• Upper limit	90°	100°	60°	100°	60°	60°
• Speed	2	2	2	2		2
• Load	6	6	6	6		6
• Extension pause	0	0	0	0		0
• Flexion pause	0	0	0	0		0
• Timer	0	0	0	0		0

### Possible values for each parameter:

	Abduction	Flexion/Extension	Rotation	Abduction + Rotation		Horizontal abduction*
• Lower limit	20° to 155°	20° to 175°	-60° to 85° extern	20° to 155°	-30° to 85° extern	-30° to 105°
• Upper limit	25° to 160°	25° to 180°	-55° to 90° extern	25° to 160°	-25° to 90° extern	-25° to 110°
• Speed	1 to 5 (from 30° to 120° per minute)					
• Load	1 to 6					
• Extension pause	0 to 900 seconds (15 minutes)					
• Flexion pause	0 to 900 seconds (15 minutes)					
• Timer	No time (00H00) to 24H00					

\* Available on Centura 5 or in option.

## USE OF THE HAND CONTROL

GB

### How to adjust the parameters of single movements



Beginning	Keys to press	Display	Remarks
To stop the unit		0 ROTATION STOP 0 60	Check if the locking switch is in the following position:
To choose the movement		30 Abd/add STOP 90 90	The display shows the new movement selected and the default settings of the upper and lower limits of this movement.
Or		0 Rotation STOP 0 60	
Or		30 FLEX/EXT STOP 89 100	
Or		30 ABD HOR STOP 50 60	The display shows "not available" if you don't have the horizontal abduction module.
To display the lower limit of the movement		30 FLEX/EXT STOP 89 100	The value blinks.
To change the lower limit if necessary		50 FLEX/EXT STOP 89 100	The new value blinks.
To validate the new value, press another key or wait more than 3 seconds	 timer speed force	50 FLEX/EXT STOP 89 100 FLEX/EXT NO TIMER FLEX/EXT SPEED 102 FLEX/EXT LOAD : ■■■■	While the value blinks press the  or  key to change if necessary.
Or to display pause		50 FLEX/EXT PAUSE HIGH 105	The pause value in upper limit of movement blinks.
To change pause value in upper limit of movement if necessary		50 FLEX/EXT PAUSE HIGH 105	The new pause value in upper limit of movement blinks.
To validate and display the pause value in lower limit of movement		50 FLEX/EXT PAUSE LOW 105	The pause value in lower limit of movement blinks.
To change the pause value in lower limit of movement if necessary		50 FLEX/EXT PAUSE LOW 105	The new pause value in lower limit of movement blinks.
To validate the new value press another key or wait more than 3 seconds. The display shows the selected movement.		30 FLEX/EXT STOP 89 100	The unit is ready to start with the new parameters.

## USE OF THE HAND CONTROL

### How to set the synchronized movement parameters



GB

Beginning	Keys to press	Display	Remarks
To stop the unit		50 FLEX/EXT STOP 15 100	Check if the locking switch is in the following position:
To select the combined movement		30 SYNC ABD ROT STOP 15 100	The indication "ABD" blinks and the display shows the values for the abduction movement. To change it, proceed as for a single movement.
To press a second time on the key		0 SYNC ABD ROT STOP 0 60	The indication "ROT" blinks and the display shows the values for the rotation movement. To change it, proceed as for a single movement.

#### Synchronization rules:

- The degrees of rotation are lower than or equal to the degrees of abduction.
- 1° of abduction means 1° of rotation.
- When the degrees of rotation are lower than the degrees of abduction, the synchronization applies to the upper degrees of the movement.

Example : abduction from 30° to 100°  
rotation from 50° to 90°



#### Comments:

- Speed, load, pauses and timer are the same for both of the movement components. The setting is the same as for a single movement.
- Pauses can be set at the lower and/or the upper limits of the abduction movement.
- You will have successive displays of abduction movement limits, or associated rotation movement, by repeatedly pressing the synchronized movement button.
- You cannot change the settings while the machine is running.

## USE OF THE HAND CONTROL










### Using Programs

The KINETEC Centura allows you to store up to 16 programs, including the type of movement, ROM, speed, load, pauses and timer.

The original parameter values of the program are empty.

These values can be modified and recorded at any time (see 'How to enter a program' p 14).

#### To select a program:

















Beginning	Keys to press	Display	Remarks
To stop the unit		30 FLEX/ext STOP 89 100	Check if the locking switch is in the following position:  or 
To access the program mode		PROGRAM 1 EMPTY	The program number blinks.
To change the program if necessary	 or 	0 PROGRAM 3 ROTATION 90	The new program number blinks.
To exit and validate the selected program		0 ROTATION STOP 0 60	The current parameters have been recorded in program 3.
To exit without validation of selected program		30 FLEX/ext STOP 89 100	Back to the starting parameters.
Start the unit		30 FLEX/ext RUN 89 100	The value change at the speed of the movement.

#### Comments:

- The values show in the 'Display' column are examples. They actually depend on the stored programs.
- The current movement parameters can be changed while using that program but no data will be stored in the original program. See the programming mode (p 14) to modify programs.

## USE OF THE HAND CONTROL

### Reading the values of a program: example SPEED




















Beginning	Keys to press	Display	Remarks
To stop the unit			Check if the locking switch is in the following position:  or 
To access the program mode			The program number blinks.
To change the program if necessary	 or 		The new program number blinks.
To read the speed value			Displaying the speed value.
After 15 seconds or after pressing on another key			
To exit and validate the selected program			The current parameters have been recorded in program 3.
Start the unit			The value change at the speed of the movement.

#### Comments:

- The values showed in the 'Display' column are examples. They actually depend on the stored programs.
- The current movement parameters can be changed while using that program but no data will be stored in the original program. See the programming mode (p 14) to modify programs.

## USE OF THE HAND CONTROL













### How to modify programs PROGRAM MODE:

Beginning	Keys to press	Display	Remarks
To switch off the unit			Check if the locking switch is in the following position: 
To press the two keys at the same time to switch the unit on	 	KINETEC CENTURA V3.0	Welcome text during 3 seconds.
Then		PROGRAM 1 EMPTY	The program number blinks.
To change the program if necessary	 or 	PROGRAM 10 EMPTY	The new program number blinks.
To choose the movement		PROGRAM 10 30 ABD/ADD 90	The display indicates the selected movement, the program number blinks again.
Or		PROGRAM 10 0 ROTATION 60	
Or		PROGRAM 10 30 FLEX/EXT 60	
Or		PROGRAM 10 30 SYNC ABD 00	
Or		PROGRAM 10 30 ABD HOR 00	The display shows "not available" if you don't have the horizontal abduction module.
To display the lower limit of the movement		PROGRAM 10 30 SYNC ABD 100	The value blinks.
To change the lower limit of the movement if necessary	 or 	PROGRAM 10 50 SYNC ABD 100	The new value blinks.
To validate the new value, press another key	 timer speed force	PROGRAM 10 50 SYNC ABD 100 PROGRAM 10 TIMER 00H15 PROGRAM 10 SPEED : 1 PROGRAM 10 LOAD : 	The value blinks. Press the  or  key to change if necessary.
Or display pauses	pause	PROGRAM 10 PAUSE HIGH 0	The pause value in upper limit of movement blinks.
To change pause value in upper limit of movement if necessary	 or 	PROGRAM 10 PAUSE HIGH 10	The new pause value in upper limit of movement blinks.

## USE OF THE HAND CONTROL

### How to modify programs PROGRAM MODE (continued)

GB

Beginning	Keys to press	Display	Remarks
To validate and display the pause lower limit of movement	pause	PROGRAM 10 PAUSE low 0	The pause value in lower limit of movement blinks.
To change the pause lower limit of movement if necessary	 or 	PROGRAM 10 PAUSE low 20	The new pause value in lower limit of movement blinks.
To validate and display of the combined rotation setting	 abd/add+rot	PROGRAM 10 0 sync rot 60	The program number blinks and the display indicates the rotation values combined with abduction.
To change the lower limit of the movement		PROGRAM 10 0 sync rot 60	The value blinks.
To change the lower limit if necessary	 or 	PROGRAM 10 10 sync rot 60	The new value blinks.
To validate and display the upper limit of the movement		PROGRAM 10 10 sync rot 60	The value blinks.
To change the upper limit if necessary	 or 	PROGRAM 10 10 sync rot 75	The new value blinks. (see page 11 for more information about combined movement)
To record the program 10	program	PROGRAM 10 Save: + clear: -	
Then		PROGRAM 10 Saving PROGRAM empty 11	The program 10 has been recorded and the display indicates the next program so you can change another program.
OR		PROGRAM 10 clearing PROGRAM empty 11	The program 10 has been cancelled and the display indicates the next program so you can change another program.
To exit program mode, switch off and switch on the unit.		KINETEC CENTURA V3.0	To use the modified program see page 12.

GB

## USE OF THE HAND CONTROL

### Comments:

- When a program has been deleted, the display shows **PROGRAM 11  
EMPTY**
- The values shown in the 'Display' column are examples. They actually depend on the stored programs.

### Program table:

Program number	Movement type	Lower limit	Upper limit	Speed	Load	Pause on upper limit	Pause on lower limit	Timer
1								
2								
3								
4								
5								
6								
7								
8								
9								
10								
11								
12								
13								
14								
15								
16								









## USE OF THE HAND CONTROL

### How to define the upper and lower movement limits

#### • At the start of a session

The MANUAL MODE is a way to set within the tolerance of a patient at the beginning of a session.

Proceed as below:

Beginning	Keys to press	Display	Remarks
Switch the unit on		<div>KINETEC CENTURA V3.0</div> <div>MOVEMENT VERIF. PLEASE WAIT</div> <div>MOVEMENT VERIF. MOTEO: M1 M2</div> <div>FLEX/EXT 50 STOP 15 100</div>	Check if the locking switch is in the following position:  or 
To select the MANUAL MODE for upper limits by continually holding pressure on the key	 continuous press	<div>FLEX/EXT</div> <div>50 MANUaL 15 100</div>	The unit is moving to the upper limit of the movement.
To set the pain level when reached, immediately press	 limit	<div>FLEX/EXT</div> <div>50 MANUaL 150 150</div>	The new upper value limit of the movement is recorded.
To select the manual mode for lower limits	 continuous press	<div>FLEX/EXT</div> <div>50 MANUaL 100 150</div>	The unit is moving to the lower limit of the movement
To set the pain level when reached, immediately press	 limit	<div>FLEX/EXT</div> <div>25 MANUaL 25 150</div>	The new lower value limit of the movement is recorded.
To start the session with the new movement limits	 start	<div>FLEX/EXT</div> <div>25 run 30 150</div>	The angle display changes with current movement.

#### Specific rules for synchronized movement:

- You can only change the upper limit of the movement through the manual mode and only beyond the synchronization point.

#### Comments:

- The values shown in the 'Display' column are examples. They actually depend on the stored programs.














GB

GB

## USE OF THE HAND CONTROL

### • During the session

The BY-PASS MODE is a way to record the pain threshold of a patient during a session.

Beginning	Keys to press	Display	Remarks
The unit is running	 start		The angle display changes with current movement. Check if the locking switch is in the following position:  or 
To select the BY-PASS MODE	 continuous press		The unit exceeds the recorded upper limit.
To set the new pain level when reached, immediately press	 limit		The new upper value limit of the movement is recorded.
To select the BY-PASS mode for lower limits	 continuous press		The unit is moving to the lower limit of the movement.
To set the new pain level when reached, immediately press	 limit		The new lower value limit of the movement is recorded.
Continue the session with the new movement limits.			The angle display changes with current movement.

### Specific rules for synchronized movement:

- You can only change the upper limit of the movement through the manual mode and only beyond the synchronization point.


### Comments:

- The values shown in the 'Display' column are examples. They actually depend on the stored programs.


## INSTRUCTIONS FOR USE

### How to use the pads

KINETEC Centura is delivered with 7 straps:  
Part number to order the complete set: 4650001397



- 4 straps on the abduction combined with the rotation or horizontal abduction splint.





- 3 straps on the abduction or flexion splint.

All these straps are used the same way (see pictures).

Do not tighten the straps too much.

To meet optimal hygiene, clean the straps after each patient.



# INSTRUCTIONS FOR USE

## Patient set up


Make sure the straps are clean.

Put the unit in the position that is the most comfortable for the patient.


GB



Position the patient in the chair in a comfortable position and supporting the affected arm.




Slide the arm supports toward the patient and put the arm in the supports. Secure pads.



Adjust the lengths:  
a - arm  
b - forearm.

## Adjusting the shoulder joint axis

- c – Vertical adjustment:
  - Loosen the 2 knobs (1).
  - With the handle (2), adjust the height of the entire mechanism.
  - Tighten the 2 knobs (1).
- d – Side to side adjustment:
  - Loosen the knob (3)
  - Slide the entire mechanism
  - Tighten the knob (3)
- e – Scapula plan choice
  - Loosen the knob (4)
  - Rotate the arm support
  - Tighten the knob (4).
- Adjust the arm rest.




## Starting the unit

- Adjustment of the rotation position:

- Press  and find the right position with the MANUAL MODE (see page 17)

- Choice of the abduction/adduction motion:


- Press  and adjust your parameters (see page 10).

- Or select a program  (see page 12).

GB

## INSTRUCTIONS FOR USE

### Adjustments for ABDUCTION/ADDUCTION with fixed ELBOW EXTENSION / FLEXION



The KINETEC Centura provides motion from 20° to 160° of abduction.




During this motion the elbow flexion settings are fixed.

#### Parts needed


• Chair

• Abduction or flexion splint


• Hand control




#### Assembling the parts




Loosen the knob (1) and slide the motor support (2) to the right or the left. Plug in the hand control.




Assemble the abduction motor (3) and tighten the screw (4). Plug in the motor.




Use the color code to assemble the forearm splint (red for right, blue for left). The assembly is secure when you hear a 'click'.



Use the color code to assemble the forearm splint (7). Tighten the knob (6).



Position arm splint (8). The assembly is secure when you hear a 'click'.



The KINETEC Centura is shown assembled for left shoulder mobilization.

# INSTRUCTIONS FOR USE


## Patient set up

Make sure the straps are clean.


Put the unit in the position that is the most comfortable for the patient.



Position the patient in the chair in a comfortable position and supporting the affected arm.




Slide the arm supports toward the patient and put the arm in the supports. Secure pads.



Adjust the lengths:  
a - arm  
b - forearm.

## Adjusting the shoulder joint axis

- c – Vertical adjustment:
  - Loosen the 2 knobs (1).
  - With the handle (2), adjust the height of the entire mechanism.
  - Tighten the 2 knobs (1).
- d – Side to side adjustment:
  - Loosen the knob (3)
  - Slide the entire mechanism
  - Tighten the knob (3)
- e – Scapula plan choice
  - Loosen the knob (4)
  - Rotate the arm support
  - Tighten the knob (4).
- Adjust the arm rest.




## Starting the unit


- Adjustment of the elbow flexion position:

- Unscrew the knob (6) and adjust the flexion as appropriate.

- Choice of the abduction/adduction motion:

- Press  and adjust your parameters (see page 10)


- Or select a program  (see page 12).



## INSTRUCTIONS FOR USE



### Adjustments for ABDUCTION/ADDUCTION with ASSOCIATED ROTATION

The KINETEC Centura provides motion from 20° to 160° of abduction associated with 120° of rotation in maximum.




#### Parts needed


- Chair
- Abduction with associated rotation splint
- Hand control




#### Assembling the parts




Loosen the knob (1) and slide the motor support (2) to the right or the left. Plug in the hand control.




Assemble the abduction motor (3) and tighten the screw (4). Plug in the motor.



Use the color code to assemble the rotation motor (red for right, blue for left). The assembly is secure when you hear a 'click'. Plug in the motor.



Use the color code to assemble the elbow splint (6). The assembly is secure when you hear a 'click'.



Use the color code to assemble the forearm splint (7). Tighten the knob (8).



The KINETEC Centura is shown assembled for left shoulder mobilization.

# INSTRUCTIONS FOR USE


## Patient set up

Make sure the straps are clean.


Put the unit in the position that is the most comfortable for the patient.



Position the patient in the chair in a comfortable position and supporting the affected arm.




Slide the arm supports toward the patient and put the arm in the supports.  
Secure pads.



Adjust the lengths:  
a - arm  
b - forearm.


## Adjusting the shoulder joint axis

- c – Vertical adjustment:
  - Loosen the 2 knobs (1).
  - With the handle (2), adjust the height of the entire mechanism.
  - Tighten the 2 knobs (1).
- d – Side to side adjustment:
  - Loosen the knob (3)
  - Slide the entire mechanism
  - Tighten the knob (3)
- e – Scapula plan choice
  - Loosen the knob (4)
  - Rotate the arm support
  - Tighten the knob (4).
- Adjust the arm rest.



## Starting the unit


- Choice of the abduction/adduction associated with rotation motion:

- Press  and adjust your parameters (see page 11).

- Or select a program  (see page 12).

# INSTRUCTIONS FOR USE

## Adjustments for EXTERNAL ROTATION motion





The KINETEC Centura provides motion from 60° of internal rotation to 90° of external rotation.


During this motion the abduction settings are fixed.

### Parts needed


- Chair
- Abduction rotation splint
- Hand control




### Assembling the parts




Loosen the knob (1) and slide the motor support (2) to the right or the left. Plug in the hand control.




Assemble the abduction motor (3) and tighten the screw (4). Plug in the motor.




Use the color code to assemble the rotation motor (red for right, blue for left). The assembly is secure when you hear a 'click'. Plug in the motor.



Use the color code to assemble the elbow splint (6). The assembly is secure when you hear a 'click'.



Use the color code to assemble the forearm splint (7). Tighten the knob (8).



The KINETEC Centura is shown assembled for left shoulder mobilization.


# INSTRUCTIONS FOR USE

## Patient set up


Make sure the straps are clean.  
Put the unit in the position that is the most comfortable for the patient.



Position the patient in the chair in a comfortable position and supporting the affected arm.




Slide the arm supports toward the patient and put the arm in the supports.  
Secure pads.



Adjust the lengths:  
a - arm  
b - forearm.

## Adjusting the shoulder joint axis

- c – Vertical adjustment:
  - Loosen the 2 knobs (1).
  - With the handle (2), adjust the height of the entire mechanism.
  - Tighten the 2 knobs (1).
- d – Side to side adjustment:
  - Loosen the knob (3)
  - Slide the entire mechanism
  - Tighten the knob (3)
- e – Scapula plan choice
  - Loosen the knob (4)
  - Rotate the arm support
  - Tighten the knob (4).
- Adjust the arm rest.




## Starting the unit

Adjustment of the abduction position:

- Press  and find the right position with the MANUAL MODE (see page 17)

• Choice of the rotation motion:

- Press  and adjust your parameters (see page 11).


- Or select a program  (see page 12).

# INSTRUCTIONS FOR USE

## Adjustments for EXTENSION / FLEXION

The KINETEC Centura provides motion from 20° to 180° of flexion.

During this motion the elbow flexion settings are fixed.





### Parts needed


• Chair

• Abduction or flexion splint


• Hand control




### Assembling the parts



Loosen the knob (1) and slide the motor support (2) to the right or the left. Plug in the hand control.




Assemble the abduction motor (3) and tighten the screw (4). Plug in the motor.




Use the color code to assemble the forearm support (red for right, blue for left). The assembly is secure when you hear a 'click'.



Use the color code to assemble the forearm splint (6). Tighten the knob (7).



Position arm splint (8). The assembly is secure when you hear a 'click'.



The KINETEC Centura is shown assembled for left shoulder mobilization.

# INSTRUCTIONS FOR USE

## Patient set up

Make sure the straps are clean.


Put the unit in the position that is the most comfortable for the patient.



Position the patient in the chair in a comfortable position and supporting the affected arm.



Slide the arm supports toward the patient and put the arm in the supports. Secure pads.




Adjust the lengths:  
a - arm  
b - forearm.

## Adjusting the shoulder joint axis


- c – Vertical adjustment:
  - Loosen the 2 knobs (1).
  - With the handle (2), adjust the height of the entire mechanism.
  - Tighten the 2 knobs (1).
- d – Side to side adjustment:
  - Loosen the knob (3)
  - Slide the entire mechanism
  - Tighten the knob (3).
- e – Shoulder depth adjustment:
  - Loosen the knob (4).
  - Slide the entire mechanism.
  - Tighten the knob (4).
- f – Scapula plan choice
  - Loosen the knob (5)
  - Rotate the arm support
  - Tighten the knob (5).


- Adjust the arm rest.



## Starting the unit

- Choice of flexion/extension:

- Press  and adjust your parameters (see page 10).

- Or select a program  (see page 12).


GB

## INSTRUCTIONS FOR USE

### Adjustments for HORIZONTAL ABDUCTION



*Available on Centura 5 or in option,  
Contact your nearest KINETEC distributor.*

The KINETEC Centura provides motion from  $-30^{\circ}$  to  $110^{\circ}$  of horizontal abduction.




#### Parts needed


- Chair
- Horizontal abduction splint
- Hand control




#### Assembling the parts



Assemble the horizontal abduction column to the unit.  
Use the colour code to turn the column in the right position.  
Fasten the screw.



Assemble the motor M1 to the column. You need to see  
the red point on the M1 part for a right side set up (blue  
for left). Several positions are possible to adjust the  
plan of the motion (horizontal to more or less inclined).



Insert the elbow support holder into the extremity of the  
M1 part. The required colour point ( red for right and  
blue for left ) needs to be positioned internally. For  
safety reason, this part cannot be removed from the  
lower part and must be so lifted when removed.



Insert the elbow splint into the elbow support holder  
respecting the colour code.

**KINETEC**


a AbilityOne

# INSTRUCTIONS FOR USE

## Patient set up

Make sure the straps are clean.

Put the unit in the position that is the most comfortable for the patient.





Position the patient in the chair in a comfortable position and supporting the affected arm.

Slide the arm supports toward the patient and put the arm in the supports.  
Secure pads.

Adjust the lengths:  
a - arm  
b - forearm.

## Adjusting the shoulder joint axis


- c – Vertical adjustment:
  - Loosen the 2 knobs (1).
  - With the handle (2), adjust the height of the entire mechanism.
  - Tighten the 2 knobs (1).
- d – Side to side adjustment:
  - Loosen the knob (3)
  - Slide the entire mechanism
  - Tighten the knob (3).
- e – Shoulder depth adjustment:
  - Loosen the knob (4).
  - Slide the entire mechanism.
  - Tighten the knob (4).
- f – Plan of the motion: from horizontal to inclined
  - Loosen the knob (5)
  - Rotate the arm support
  - Tighten the knob (5).



- Adjust the arm rest.

## Starting the unit


- Choice of Horizontal abduction:

- Press  and adjust your parameters (see page 10).

- Or select a program  program (see page 12).

## OPTIONS


GB



Horizontal abduction Module



Seat height adjuster



Shoulder fixation



Head rest



Recharger



Scales kit



Transport bag

## PRODUCT INFORMATION

### MAINTENANCE

After 2,000 hours of operation, KINETEC Centura requires a few lubrication and maintenance operations (lubrication of the joints, pointer stops and ballscrews). The need for maintenance is indicated by display of the message SERV. MOTOR when the system is switched on.

Despite that warning, you can continue to use your KINETEC by pressing START, but you should contact your nearest KINETEC technician to have the maintenance operations conducted as soon as possible.

When the system is no longer in operating order, kindly return it to us with the accessories for destruction.

### TROUBLE SHOOTING

A spare parts list and technical catalog are available to you on request from your KINETEC distributor.

If, after connecting the power supply cable to the power supply and switching on KINETEC Centura:

- The display does not indicate any information:
  - Check that the electrical socket is live using another device.
  - Replace the fuse(s) of the connector with fuses of the same type and caliber: 2 fuses T 750 mA 250V (6.3 x 32) (KINETEC order: 4610007434).
  - If the display still does not indicate any information, contact your nearest KINETEC technician.

If, after switching on your KINETEC:

- Your KINETEC does not work and the display indicates 50 STOP 25 115, Press START again.

- Your KINETEC still does not function:

Contact your nearest KINETEC technician.

- Your KINETEC does not function and the display indicates:

ANGULAR POSI.: angle measurement function failure,  
or NO MOVEMENT: no movement,  
or BAD WAY: motor rotation failure,  
or LOAD MAXI: abnormal consumption,  
or POWER SUPPLY: power failure;

Contact your nearest KINETEC technician if the same message is displayed after having switched the device off, then on, and started it by pressing START.

### CLEANING

Before conducting any cleaning operation, SWITCH OFF the unit and disconnect the power supply.


Use a DISINFECTANT (PROPANOL/ISOPROPANOL or ALDEHYDE-based solution). Spray the disinfectant on the SURFACES (plastic shells and metal components).

In order to ensure optimal hygiene, you are advised to clean the covers for each new patient. All the consumables enable hazard-free disposal.

#### Recommendations to obtain a maximum hygiene of the pads:

- Sterilization of the pads ( if necessary) :  
Sterilizer at 134 °C during 18 minutes.
- Desinfection of the pads:  
Washing at 30 ° C with use of a disinfection solution during the rising cycle.  
Example of product which can be used :  
Solution " Baclinge" at 0.125 % or "Souplanios" at 0.125% from ANIOS Laboratory.  
A complete list of distributors in your country is available on request.

GB



## PRODUCT INFORMATION

### TECHNICAL SPECIFICATIONS

GB

Product  
Weight: 22 Kg /48lb  
Splint dimensions:  
56x100x76cm / 22"x39"x30"  
Angular limits: see page 2  
Speeds: from 30 to 120°/min  
Patient height: from 1,40 to 2m  
47" to 6'7"

Electricity  
Power supply: 100-240 V~  
Frequency: 50-60 Hz  
Power consumption: 50 VA  
Device of type B class I  
IP 20.  
Fuse T 750mA 250V 6.3x32mm  
KINETEC order: 4610007434

Environment  
- Storage/transport conditions:  
Temperature: -40 to 70°C / -40 to 160°F  
Relative humidity: up to 90%  
- Operating conditions:  
Room temperature: 10 to 40°C / 50 to 105°F  
Relative humidity: up to 80%

### SYMBOLS USED

	TYPE B device (protection against electric shocks)
	Caution (consult the accompanying documents)
	STOP (power off)
	ON (power on)
	Start movement
	Stop movement
	Program access
	Speed
	Timer
	Force
	Pause
	Increase
	Decrease

	Lower limit
	Upper limit
	Flexion movement
	Rotation movement
	Abduction movement
	Combined movement
	Horizontal abduction movement
	Hand control locked
	Hand control unlocked
	Hand control half locked
	Switch on LED and defect signal when the LED blinks
	Alternative current

### WARRANTY

The KINETEC warranty is strictly limited to the replacement free of charge or repair in the plant of the component or components found to be defective.

KINETEC guarantees its joint passive mobilization systems for 1 year against all defects of manufacture from the date of purchase by the consumer.

KINETEC is the only organization able to assess the application of the warranty to its systems.

The warranty will be considered null and void if the device has been used abnormally or under conditions of use other than those indicated in the user's manual.

The warranty will also be considered null and void in the event of deterioration or an accident due to negligence, inappropriate surveillance or inappropriate maintenance, or due to transformation of the equipment or an attempt to repair the equipment.

**KINETEC**

AbilityOne

# DECLARATION "CE" DE CONFORMITE "EC" DECLARATION OF CONFORMITY

Nous  
We

**AbilityOne Kinetec**

(nom du fabricant ou de son mandataire dans la Communauté Européenne)  
(Name of manufacturer or his representative established within the European Community)  
Zone Industrielle de Tournes-Cliron - TOURNES - F 08090

déclarons sous notre seule responsabilité que le produit  
declare on our own responsibility that product

## KINETEC d'épaule CENTURA SHOULDER KINETEC model CENTURA

(nom, type ou modèle, (éventuellement n° de lot, d'échantillon ou de série, sources et nombre d'exemplaires))  
(Brand name, model, (N° of lot or series if necessary))

satisfait aux dispositions des Directives du Conseil :  
complies with the assessment criteria of Council's Directives:

- n° 93/42/CEE du 14 juin 1993 - DISPOSITIFS MEDICAUX (selon l'annexe II)  
n° 93/42/CEE of June 14, 1993 - MEDICAL DEVICES (Annex II)  
obligatoire à partir du 14 Juin 1998 / mandatory date : June 14, 1998
- n° 89/336/CEE du 3 mai 1989 - COMPATIBILITE ELECTROMAGNETIQUE  
n° 89/336/CEE of May 3, 1989 - ELECTROMAGNETIC COMPATIBILITY  
obligatoire à partir du 1er janvier 1996 / mandatory date : January 1, 1996

et que le système qualité du fabricant est conforme aux normes ou autres documents normatifs suivants :  
and the quality system of the manufacturer is in conformity with the following standard (s) :

EN 46001/96 - ISO 9001/94  
Certificat n° 0214/46001/9001/1 délivré par le Gmed le 15/10/97  
Certificate n° 0214/46001/9001 established by Gmed on 15/10/97

Attestation de conformité à l'annexe II.3 délivrée par le Gmed N°0214/B2P3/1  
Certificate for approval of full Quality Assurance System - Annexe II.3 - N° 0214/B2P3/1 established by Gmed

Information complémentaire : Fiche Produit Centura du 04/04/2002  
Additional information : Product file Centura dated 04/04/2002

Année d'apposition du Marquage CE : 2000  
Year in which CE mark was affixed : 2000

Tournes, le 4 Avril 2002  
Tournes, April 4<sup>th</sup>, 2002  
(lieu et date / place and date)




**P. SONNET**  
Directeur Général / General Manager  
(Nom, titre et signature du signataire autorisé)  
(Name, function and sign of authorised people)

CENTURA - REV2


## SELECT YOUR LANGUAGE

GB

Beginning	Keys to press	Display	Remarks
Switch ON the unit		<div>KINETEC CENTURA V3.0</div> <div>VERIFICATION MVT ATTENDEZ S.V.P</div> <div>VERIFICATION MVT MOTEUR: M1 M2</div> <div>FLEX/EXT 30 STOP 89 100</div>	<p>Check if the locking switch is in the following position</p> 
Press the 2 keys in the same time	speed force	Language FRENCH	The display indicate the language selectionned
To change the language	+ or -	Language ENGLISH	<p>The English language is selectionned.</p> <p>English French German Italian Spanish</p>
To validate the new language.	limit	OK Switch on/off	To exit and confirm the new language, switch OFF and Switch ON the unit..

## Pads for Kinetec® Centura® Shoulder CPM Machine:

5315-0130




### COMPONENTS


The following pads and straps (labeled at left) are included in the pad kit:

- a. elbow pad
- b. wrist pad
- c. hand strap
- d. biceps pad


### • Components



- Abduction/adduction
- Rotation
- Synchronized abduction/adduction with rotation



### • Flexion/extension



### • Hand strap and wrist pad attachment

### SETTING UP THE MACHINE

Each set-up requires only some of the pads, as follows:

- Set-up for abduction/adduction requires pads **a**, **b**, and **c**.
- Set-up for rotation requires **a**, **b**, and **c**.
- Set-up for synchronized abduction/adduction with rotation requires **a**, **b**, and **c**.
- Set-up for flexion/extension requires **b**, **c**, and **d**.

Attach the pads and straps to the CPM machine as follows:

- a. Elbow pad**—Position it in the elbow support with the seam distal and the straps and D-rings toward the elbow support. Thread the straps and D-rings through the slots in the elbow support.
- b. Wrist pad**—Position it on the wrist support with the strap and D-ring toward the wrist support. Thread the strap and D-ring through the slots in the wrist support.
- c. Hand strap**—Turn the hand support so the rod is up. With the hook material on the strap facing up, thread the D-ring end of the strap under the rod.
- d. Biceps pad**—Position it on the biceps support with the strap and D-ring toward the biceps support. Thread the strap and D-ring through the slots in the biceps support.

After patient has been placed in the supports, secure straps by threading them through the D-rings and securing the hook closures. **Be sure the straps are snug but not so tight that they impair circulation.**

**The pads are intended for single-patient use only.**



Certificate 2/3 N° : 0214 / B2P3 / 1

Delivered at Fontenay-aux-Roses on November 12<sup>th</sup>, 1997


## EC CERTIFICATE

APPROVAL OF FULL QUALITY ASSURANCE SYSTEM  
ANNEX II point 3 DIRECTIVE 93/42/EEC concerning medical devices

Device(s) category :	Immobilization, mobilization and testing devices of the human joints
Identification of device(s) :	See EC declaration of conformity of the manufacturer authenticated by G-MED
Manufacturer(name and address) :	SMITH & NEPHEW - KINETEC S.A. Zone Industrielle de Tournes - Cliron BP 19 F - 08090 TOURNES
E.U. Responsible :	SMITH & NEPHEW - KINETEC S.A. Zone Industrielle de Tournes - Cliron BP 19 F - 08090 TOURNES

The G-MED certifies that, on the basis of the results contained in the file referenced 30273304, the quality assurance system - for design, production and final inspection - of medical devices designated above is in conformance with the requirements of the directive 93/42/EEC, annex II point 3.

This certificate is valid until : November 12<sup>th</sup>, 2002 (included)



E. GRIMAUD  
G-MED General Manager

In compliance with clause 17, CE marking shall be accompanied by our identification number 0459  
This certificate is issued under the following conditions

- 1 It is valid only to the designated device(s) listed above
- 2 Copies are available upon request by the manufacturer

REVISION: A/ AAT2DAN2/LE 1994-12-21

33 AVENUE DU GÉNÉRAL LECLEF - BP 19 - 08090 TOURNES  
TEL INTER : +33 1 40 95 53 54 - FAX INTER : +33 1 40 95 53 55  
Groupeement d'entreprises pour l'évaluation des dispositifs médicaux