ARTROMOT[®]-K1



Service Manual

Table of contents

4 11:-4-....

I. HISTORY	2
2. Purpose	2
3. General	2
4. Packing and unpacking the ARTROMOT®-K1	5
5. Block diagram of the electronic parts ARTROMOT®-K1	6
6. Bill of material for service parts ARTROMOT®-K1	7
7. Figure for bill of material	8
8. Special function Service menu for ARTROMOT®-K1Comfort and Standard version	10
9. Special function Service menu for ARTROMOT®-K1Classic	12
10. How to perform repairs	14
11. Checklist of safety and function test ARTROMOT®-K1Comfort and Standard version	15
12. Checklist of safety and function test ARTROMOT®-K1Classic	16

1. History

Revision	Date	Name	Change
1	10.02.2006	S. Herr	Service Manual created
2	14.07.2006	S. Herr	Pos. 25, chapter 10
3	14.08.2006	S.Herr	Pos. 20, Chapter 9, 12

2. Purpose

The purpose of this Service Manual is to help you make simple repairs on the ARTROMOT®-K1. Only authorized staff may perform repairs and maintenance as the manufacturer's warranty and liability would otherwise be invalidated. Only original parts may be used for servicing in accordance with the attached spare parts list.

3. General

3.1 Electronics, connection cables

No plugs may be connected or disconnected while the unit is switched on. Always switch the ARTROMOT®-K1 off before connecting or disconnecting a plug.

The locks for spiral cable for the hand-held programming unit have to be closed at all times.

When you assembling with electronic parts make sure to use ESD (Electro Static Discharge) equipment.

Only original chip cards may be used. Insert the chip cards so that ARTROMOT@ is visible.

If you have to exchange any of the printed circuit boards including the knee electronics, hand held programming unit, motor control or power supply electronics you have to perform a calibration.

Possible errors: Errors will be displayed on the handheld programming unit as follows shown:



ERROR XX

(XX = Number of the error)

1 Potentiometer error:

Wrong angle information

- -> Check the femur settings
- -> Replace knee electronics (Pos. 5)
- -> Replace motor control (Pos. 25)

2 Failure at the potentiometer:

Connection to the potentiometer is interrupted

- -> Replace the spiral cable of the potentiometer (Pos. 24)
- -> Replace knee electronics (Pos. 5)
- -> Replace motor control (Pos. 25)

3 Motor driver error:

The motor driver IC reported an error

-> Replace motor control (Pos. 25)

4 Motor error:

The motor did not turn properly

- -> Replace motor control (Pos. 25)
- -> Replace the motor (Pos. 26)

5 Motor over current:

The current for the motor exceeded the maximum limit

- -> Check the mechanics (Pos. 30)
- -> Replace motor control (Pos. 25)
- -> Replace the motor (Pos. 26)

6 Motor control error:

Internal error in the motor control.

-> Replace motor control (Pos. 25)

7 Eprom access error:

Memory error in the access of the EPROM

-> Replace the hand held programming unit (Pos. 20)

8 CPM ROM error:

Memory error in the motor control

-> Replace motor control (Pos. 25)

9 Communication:

Communication to the motor control is not possible

- -> Check spiral cable and connector
- -> Replace the hand held programming unit (Pos. 20)
- -> Replace motor control (Pos. 25)

10 General error in the motor control:

Unknown error in the motor control

-> Replace motor control (Pos. 25)

11 Motor enable timeout

Motor could not be enabled in time

-> Replace motor control (Pos. 25)

12 Invalid parameter motor error:

Motor has received a invalid parameter

- -> Replace motor control (Pos. 25)
- -> Replace the hand held programming unit (Pos. 20)

13 Stop release error:

The motor could not be released

- -> Replace motor control (Pos. 25)
- -> Replace the hand held programming unit (Pos. 20)

14 Unexpected motor Stop:

- -> Check spiral cable and connector
- -> Replace motor control (Pos. 25)

15 Motor disabled:

Motor control disabled the motor.

-> Replace motor control (Pos. 25)

16 Wrong command in the motor:

- -> Replace motor control (Pos. 25)
- -> Replace the hand held programming unit (Pos. 20)

17 5V supply error:

5V supply of motor control not sufficient

-> Replace motor control (Pos. 25)

18 Initialise error real time clock:

-> Replace the hand held programming unit (Pos. 20)

19 Communication error real time clock:

-> Replace the hand held programming unit (Pos. 20)

20 Error real time clock:

-> Replace the hand held programming unit (Pos. 20)

21 Range exceeded:

The measured angle is out of the range of motion

-> Replace motor control (Pos. 25)

22 ROM error in the hand held programming unit:

Memory error in the hand held programming unit

-> Replace the hand held programming unit (Pos. 20)

23 Invalid parameter:

Internal error in the hand held programming unit

-> Replace the hand held programming unit (Pos. 20)

24 24V supply error motor control:

Error in the 24V supply in the motor control

- -> Replace motor control (Pos. 25)
- -> Replace the power supply electronics (Pos. 28)

25 Bus error:

Bus system error

- Replace the spiral cable of the hand held programming unit
- -> Replace the hand held programming unit (Pos. 20)
- -> Replace motor control (Pos. 25)

26 24V supply hand held programming unit:

24V supply of the hand held programming unit is defective

-> Replace the hand held programming unit (Pos. 20)

27 5V supply hand held programming unit:

5V supply of the hand held programming unit is defective

-> Replace the hand held programming unit (Pos. 20)

28 3.3V supply hand held programming unit:

3.3V supply of the hand held programming unit is defective

-> Replace the hand held programming unit (Pos. 20)

29 Calibration:

The calibration data in the motor control are wrong.

-> Perform a calibration (see chapter 8.1)

30 Calibration error:

- -> Repeat the calibration (see chapter 8.1)
- -> Replace knee electronics (Pos. 5)
- -> Replace motor control (Pos. 25)

31 Calibration timeout:

-> Replace motor control (Pos. 25)

32 Motor enable error:

The motor could not be enabled

-> Replace motor control (Pos. 25)

33 Motor disable error:

The motor could not be disabled

-> Replace motor control (Pos. 25)

34 Motor stop error:

Motor stop command timeout error

-> Replace motor control (Pos. 25)

35 Configuration error:

Invalid configuration of the hand held programming unit -> Replace the hand held programming unit (Pos. 20)

36 Parameter not valid chipcard

-> Replace the chipcard (Pos. 21)

37 Checksum error chipcard:

-> Replace the chipcard (Pos. 21)

38 Unknow chipcard error:

-> Replace the hand held programming unit (Pos. 20)

39 A chipcard of another ARTROMOT® product got

- -> Format the chipcard; Press START
- -> Replace the chipcard (Pos. 21)

40 Unkown chipcard:

-> Replace the chipcard (Pos. 21)

41 Chipcard is missing:

No chipcard is insert

- -> Insert a chipcard
- -> Replace the chipcard (Pos. 21)

42 Chipcard write error:

The chipcard is not writable

- -> Replace the chipcard (Pos. 21)
- -> Replace the hand held programming unit (Pos. 20)

43 Chipcard read error

Chipcard is not readable.

- -> Replace the chipcard (Pos. 21)
- -> Replace the hand held programming unit (Pos. 20)

44 Chipcard verify

- -> Replace the chipcard (Pos. 21)
- -> Replace the hand held programming unit (Pos. 20)

45 Wrong product combination:

Mixup between non compatible device and hand held programming unit

-> Use the correct hand held programming unit (Pos. 20)

46 Handset error internal communication:

Invalid interchip communication inside the hand held programming unit

-> Replace the hand held programming unit (Pos. 20)

47 Internal communication error motor control:

Internal communication error motor control.

-> Replace motor control (Pos. 25)

48 User stoped the special function

49 Unknown error in the motor control:

-> Replace motor control (Pos. 25).



- -> Insert a chipcard
- -> Replace the chipcard (Pos. 21)



Chipcard error

- -> Replace the chipcard (Pos. 21)
- -> Replace the hand held programming unit (Pos. 20)

3.2 Mechanics

The movable screws should not be completely unscrewed when adjustments are being made. Make sure that the movable screws are tightened for operation and transport.

The frame is unstable:

Possible cause: Bolt / screws missing or loose. Tighten the screws / bolts.

3.3 Others

Do not clean the housing or the support with grease or oil.

No solvents may be used when cleaning the ARTROMOT®-K1.

4. Packing and unpacking the ARTROMOT®-K1

The following settings must be made to transport the ARTROMOT®-K1:

Set the packing setting in the menu or move the device in a position of EXTENSION 0 degrees.

Switch off the ARTROMOT®-K1.

Remove the power cord and disconnect the hand held programming unit.

Only use original packaging for transport.

Put the hand-held programming unit into the specified cutout of the styrofoam and fix it with adhesive tape.


Set the femur lenght on maximum and the lower leg setting on 45 cm.

Set the angle joint horizontal.

Move the two styrofoam parts on the device First put the power cord on the bottom of the box and then the ARTROMOT®-K1 with the two styrofoam parts.


Hand-held programming unit

Power cord



5. Block diagramm of the electronic parts ARTROMOT-K1

Hand-held programming unit





6. Bill of material for service parts ARTROMOT®-K1

Position	Description	Ordernumber
1	Tigh support complete with belt	2.0013.073
	Tigh support	2.0013.294
	Belt with belt strap	2.0013.059
2	Clamping lever	2.0037.203
	Pin	2.0013.352
	Circlip	DIN6799D2,3A2
4	Fixation screw	GN534-32-M5sw
	Washer	DIN125D6A2
	Countersunk screw	DIN933M5x20A2
	Distance disk	DIN988D6x12x0,3
5	Knee electronics	2.0037.023
6	Lower leg support complete with belts	2.0034.074
	Lower leg support	2.0013.291
	Belt with belt strap	2.0013.059
7	Clamping lever	2.0037.202
	Pin	2.0013.352
	Circlip	DIN6799D2,3A2
8	Fixation screw	GN534-32-M5sw
	Washer	DIN125D6A2
	Countersunk screw	DIN933M5x20A2
	Distance disk	DIN988D6x12x0,3
10	Belt loop	2.0013.087
11	Footrest complete with plate, clamping lever, belt	2.0037.032
	Footrest	2.0013.300
	Footplate	2.0037.168
	Oval plate	2.0013.298
12	Fixation screw	GN534-40-M6sw
	Rubber puffer complete for angle joint	0.0037.202
13	Clamping lever	2.0037.204
15	Connection hand held programming unit	2.0037.004
16	Power switch (ON/OFF) with connection	0.0034.245
17	Fuses cap	0.0034.246
	Fuses 1 AT	0.0000.005
20	Hand held programming unit with spiral cable	
	ARTROMOT®-K1 Standard	0.0037.050
	Hand held programming unit with spiral cable	
	ARTROMOT®-K1 Standard with chipcard	0.0037.051
	Hand held programming unit with spiral cable	
	ARTROMOT®-K1 Comfort	0.0037.054
	Hand held programming unit with spiral cable ARTROMOT®-K1 Comfort with chipcard	0.0037.055
	Hand held programming unit with spiral cable ARTROMOT®-K1 Classic	0.0037.061
	Protection for hand held programming unit	0.0037.103
	7	


21	ARTROMOT® -K1 chipcard	0.0037.035
22	Housing cover K1 with sticker	0.0037.200
23	Housing cover small	0.0037.201
24	Spiral cable of the Potentiometer	0.0037.007
25	Motor control	2.0037.901RevB
26	Motor	0.0037.010
27	Wire set	2.0037.005
28	Power supply electronics	0.0034.244
29	Lip	2.0037.105
	Wire for lip	0.0031.300
30	Spindel complete	0.0037.203
	Power cord EU version	0.0034.118
	Power cord USA version	0.0034.011
	Option: Frame adapter add-on-kit	0.0037.204

7. Figure for bill of material





Position 30



8. Special function Service Menu → C ARTROMOT®- K1 Comfort and Standard version

Function of service Menu

Calibration



Display contrast



Error log



Device runtime



Entering the service menu:

Press the menu key until Service Menu shows up.



ARTRMOT®-K1 Standard Version menu level 3. ARTRMOT®-K1 Comfort Version menu level 5.

Press



for 5 seconds.



is flashing on the display.

The display will change and show: Entering code.



For the code press: 1 3 2 4

Now you see the symbols of the service menu









8.1 Calibration



ATTENTION!

Before you do a calibration switch the device OFF and ON.

Femur setting: 49 cm.

Press the symbol calibration




Display: ← ■ ■ ■

Press + or - to move the ARTROMOT®-K1 to 0 degrees.

Press START, the calibration starts automatically. The device will reach both maximum points and move between -10 bis 120 degrees with different speed.

Wait until the ARTROMOT®-K1 stops.

If the calibration was succesful the device stops at 0 degrees and show following symbols on the display:





Press STOP twice to leave the service menu. Finally, a safety and fuction test has to be performed (see chapter 11).

8.2 Display contrast



Press the symbol display contrast







Press + or - to set up the requested display contrast. You can set the display contrast from 0 – 100%.

Press STOP twice to save the settings and leave the service menu.

8.3 Error log

Press the symbol error log

You will find following information on the display:

Upper line: Number of the current showed error message and the total number of the saved error messages.

Lower line: Error message

Left side: The symbol of the causer.



= Hand held programming unit



= Motor

Press + or - to see the entries of the error log. Press STOP twice to leave the service menu.

General note to the error log:

Entries are always in english.

The entries are ordered by causer and not in temporal order.


8.4 Device runtime



Press the symbol device runtime



The display shows the device runtime



Press STOP twice to leave the service menu.

9. Special function Service Menu → C ARTROMOT®- K1 Classic

Function of service Menu

Calibration



Display contrast



Error log



Device runtime



Entering the service menu:

Switch off the ARTROMOT®-K1.

When you switch on the device press simultaneously

the Extension control



and the Flexion control



Display will show following symbol: Entering code.



For the code press the control as shown below:









Now you see the symbols of the service menu:









9.1 Calibration



ATTENTION!

Before you do a calibration switch the device OFF and ON.

Femur setting: 49 cm.

Press the Extension control.



Display:





Keep on pressing the control

Extension



(ARTROMOT®-K1 move in direction Extension)

Flexion



(ARTROMOT®-K1 move in direction Flexion)

until the ARTROMOT®-K1 reach 0 degrees.

Press START, the calibration starts automatically. The device will reach both maximum points and move between -10 bis 120 degrees with different speed.

Wait until the ARTROMOT®-K1 stops.

If the calibration was succesful the device stops at 0 degrees and show following symbols on the display:

Display:





Press STOP twice to leave the service menu.

Finally, a safety and fuction test has to be performed (see chapter 12).

9.2 Display contrast



Press the speed control.



Display:



XX %

Press the control

Extension



(value decrease)

Flexion



(value increase)

to set up the requested display contrast. You can set the display contrast from 0 - 100 %

Press STOP twice to save the settings and leave the service menu.

9.3 Error log Fcc

Press the pause control

Z

You will find following information on the display:

Upper line: Number of the current showed error message and the total number of the saved error messages.

Lower line: Error message

Left side: The symbol of the causer.



= Hand held programming unit



= Motor

Press the control

Extension



(last entry)

Flexion



(next entry)

to see the entries of the error log.

Press STOP twice to leave the service menu.

General note to the error log

Entries are always in english.

The entries are ordered by causer and not in temporal order.

9.4 Device runtime



Keep on pressing the Flexion control.



The display shows the device runtime

Display:



XX (XX = Runtime in hours).

Press STOP twice to leave the service menu.

10. How to perform repairs

10.1 How to remove the housing cover (Pos. 22 + 23).

Move ARTROMOT®-K1 in a position approximately 80 degrees.

Turn of the power OFF at the ARTROMOT®-K1 and remove the power cord.

Move the ARTROMOT®-K1 in a stable side position.

Loosen the 4 outside torx screws to remove the housing cover K1 with sticker (Pos. 22).

Loosen the 2 torx screws to remove the housing cover small (Pos. 23).

If you have to exchange any of the printed circuit boards including the knee electronics, hand held programming unit, motor control or power supply electronics you have to perform a calibration.

10.2 How to exchange the motor control (Pos. 25).

ATTENTION!

When you assembling with electronic parts make sure to use ESD (Electro Static Discharge) equipment.

Remove the housing cover, see chapter 10.1.

Pull out the connectors of the motor control.

Loosen the four screws.

Exchange the defective motor control and fix it with the screws.

Put back in the connectors in the same position.

Rebuild the housing cover.

Perform a calibration.

Finally, a safety and function test has to be performed.

10.3 How to exchange the power supply electronics (Pos. 28).

ATTENTION!

When you assembling with electronic parts make sure to use ESD (Electro Static Discharge) equipment.

Remove the housing cover, see chapter 10.1.

Pull out the connectors of the power supply electronics.

Loosen the four screws.

Exchange the defective electronics and fix it with the screws.

Put back in the connectors in the same position.

Rebuild the housing cover.

Perform a calibration.

Finally, a safety and function test has to be performed.

11. Checklist of safety and function test ARTROMOT®- K1 Comfort and Standard version

Safety test		Measured value	Date/ Signature
Protective earth conductor resistance	≤ 0,1 Ohm	Ohm	
Ground leakkage current EN 60601 / IEC 601/ VDE 0751	≤ 500 µA	μА	
Or			
Ground leakkage current as in UL 2601	≤ 300 µA	μΑ	

Check the angle in position 0 degrees. Tolerance +/- 5 degrees. Check the angle in position 60 degrees. Tolerance +/- 5 degrees. Check the angle in position 100 degrees. Tolerance +/- 5 degrees. 3. Check the emergency-off function. Start the ARTROMOT® K1 in any mode. Press any key, the ARTROMOT® K1 will stop immediately. Check this for all keys. 4. Set up the special function "new patient" →0← press to activate it. Press START. The ARTROMOT® -K1 will move automatically to 30 degrees and stop there.	Error
Check the angle in position 60 degrees. Tolerance +/- 5 degrees. Check the angle in position 100 degrees. Tolerance +/- 5 degrees. 3. Check the emergency-off function. Start the ARTROMOT® K1 in any mode. Press any key, the ARTROMOT® K1 will stop immediately. Check this for all keys. 4. Set up the special function "new patient" →0← press to activate it. Press START. The ARTROMOT® -K1 will move automatically to 30 degrees and stop there. 5. Check the set values of "new patient" →0← Extension -> Display: 25° Flexion -> Display: 35° Speed -> Display: 50 % Reverse -> Display: 25	
Start the ARTROMOT® K1 in any mode. Press any key, the ARTROMOT® K1 will stop immediately. Check this for all keys. 4. Set up the special function "new patient" →0← press to activate it. Press START. The ARTROMOT® -K1 will move automatically to 30 degrees and stop there. 5. Check the set values of "new patient" →0← Extension -> Display: 25° Flexion -> Display: 35° Speed -> Display: 50 % Reverse -> Display: 25	
4. Set up the special function "new patient" →0← press to activate it. Press START. The ARTROMOT® -K1 will move automatically to 30 degrees and stop there. 5. Check the set values of "new patient" →0← Extension -> Display: 25° Flexion -> Display: 35° Speed -> Display: 50 % Reverse -> Display: 25	
Extension -> Display: 25° Flexion -> Display: 35° Speed -> Display: 50 % Reverse -> Display: 25	
Flexion -> Display: 35° Speed -> Display: 50 % Reverse -> Display: 25	
Speed -> Display: 50 % Reverse -> Display: 25	
Reverse -> Display: 25	
Therapy Duration -> Display: 00:00	
6. Start the ARTROMOT® -K1 in the motion range between –10 to 120 degrees.	
Set the speed = to 100%.	
Both extreme points should be reached within 45 – 65 seconds.	
7. All special funtion are deactivated. Select a therapy time of 3 minutes.	
At the end of the 3 minutes the ARTROMOT®-K1 will stop automatically in the middle position.	

12. Checklist of safety and function test ARTROMOT®- K1 Classic

Safety test		Measured value	Date/ Signature
Protective earth conductor resistance	≤ 0,1 Ohm	Ohm	
Ground leakkage current EN 60601 / IEC 601/ VDE 0751	≤ 500 µA	μA	
Or			
Ground leakkage current as in UL 2601	≤ 300 µA	μА	

Function test	OK	Error
Switch on the ARTROMOT® -K1 Classic. Press the two outer control (Extension/ Flexion) simultaneously.		
Display: Software version VX.X XX.XXX.XX (X = optional) Keep on pressing.		
Display: ARTROMOT K1 Classic		
2. The maximum range of motion for Extension/ Flexion is -10 to 120 degrees.		
Check the angle in position 0 degrees. Tolerance +/- 5 degrees.		
Check the angle in position 60 degrees. Tolerance +/- 5 degrees.		
Check the angle in position 100 degrees. Tolerance +/- 5 degrees.		
3. Check the emergency-off function.		
Start the ARTROMOT® K1 Classic in any mode.		
Press any control, the ARTROMOT® K1 will stop immediately.		
Check this for all controls.		
4. Checking the "PAUSE" function.		
Select the following settings:		
Extension = 10 degrees		
Flexion = 90 degrees		
Speed = 50%		
Pause = 10 seconds		
Start the ARTROMOT®-K1 Classic.		
At the reversal points (at 10 / 90 degrees) the pause must be lasting 10 seconds		
(+/- 2 seconds).		
The speed is substantially lower than in a 100% setting.		
5. Start the ARTROMOT® -K1 Classic in the motion range between –10 to 120 degrees.		
Set the speed to 100%.		
Both extreme points should be reached within 45 – 65 seconds.		



