



Comfort and Support
You Can Count On



L0637 CODE APPROVED

Like the Evergreen™ LSO LoPro, the Evergreen™ LSO was designed to treat a variety of indication in the lumbar spine. Supportive side panels offer additional coverage if desired.

FEATURES & BENEFITS



Independent Compression

With a 4 to 1 mechanical advantage, independent drawstrings tighten the upper and lower sections on the lower back offering direct compression for the relief of low back pain.




Materials

Breathable mesh fabric adds comfort for long-term wear. Evergreen™ products may be worn next to the skin, or over clothing.



Comfortable Tabs

The large pull tabs are easy to locate, grab, and pull.



Sizing Info

The Evergreen™ Line fits patients with waists between 21 and 57 inches.

PRODUCT INFORMATION

SIZE	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE	XX-LARGE
REF	36-3702	36-3703	36-3704	36-3705	36-3706	36-3707
RANGE	21-27 in. • 53-69 cm	26-32 in. • 66-81 cm	31-37 in. • 79-94 cm	36-42 in. • 91-107 cm	41-47 in. • 104-119 cm	46-57 in. • 117-145 cm

XX-SMALL • 36-3701 • 15-21 in. (38-53 cm)