



Comfort and Support
You Can Count On



L0456 CODE APPROVED

The Evergreen™ TLSO provides motion restriction in the thoracic and lumbar spine, treating a wide variety of indications. The back panel and the shoulder straps provide the stability needed to create a dynamic environment for healing patients with complex spinal conditions.

FEATURES & BENEFITS



Independent Compression

With a 4 to 1 mechanical advantage, independent drawstrings tighten the upper and lower sections on the lower back offering direct compression for the relief of low back pain.




Materials

Breathable mesh fabric adds comfort for long-term wear. Evergreen™ products may be worn next to the skin, or over clothing.



Comfortable Pull Tabs

The large pull tabs are easy to locate, grab, and pull.



Sizing Info

The Evergreen™ Line fits patients with waists between 21 and 57 inches.

PRODUCT INFORMATION

SIZE	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE	XX-LARGE
REF	36-5602	36-5603	36-5604	36-5605	36-5606	36-5607
RANGE	21-27 in. • 53-69 cm	26-32 in. • 66-81 cm	31-37 in. • 79-94 cm	36-42 in. • 91-107 cm	41-47 in. • 104-119 cm	46-57 in. • 117-145 cm

XX-SMALL • 36-5601 • 15-21 in. (38-53 cm)