





Product Comparison

the medcom group, ltd

 Click product names to view! 	ThermX	Game Ready	Vascutherm 5	Nice 1	Breg Polar Care Wave	Breg Polar Cube	Cryo Cuff (Gravity)	Cryo Cuff (IC)	Dr. Active CCT	DonJoy Iceman	ARTEK Sport
Rental Option	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✗
RX Needed	✓	✓	✓	✓	✓	✓	✗	✓	✗	✓	✓
Iceless	✓	✗	✓	✓	✗	✗	✗	✗	✓	✗	✗
Cold Temperature Range	38°F - 55°F	35°F - 50°F	43° - 49°F	43°F - 59°F	45°F-50°F	35°F - 50°F	35°F - 50°F	35°F - 50°F	N/A	35°F - 50°F	35°F - 50°F
Wraps Available	Knee, Hip, Shoulder, Lumbar, Universal, Ankle	Knee, Articulated Knee, Shoulder, Hip, Ankle	Knee, Shoulder, Hip, Elbow, Foot/Ankle, Hand, Lumbar, Cervical	Knee, Articulated Knee, Shoulder, Hip, Ankle, Lumbar, Butterfly	Knee, Shoulder, Universal	Knee, Shoulder, Hip, Ankle, Butterfly	Knee, Ankle, Elbow	Knee, Ankle, Elbow	Knee, Shoulder, Ankle, Wrist	Shoulder, Universal	Universal
Weight (empty)	14 lbs	8 lbs	18 lbs	10 lbs	5 lbs	8 lbs	5 lbs	5 lbs	2 lbs	5 lbs	4 lbs
Dimensions	10" x 9" x 15"	5.5" x 16" x 8.5"	6.6" x 12" x 13.2"	8" x 7.64" x 8"	12" x 12"	12" x 12"	N/A	N/A	N/A	12.75" x 12.5" x 12.25"	10" x 8" x 8"
Compression Option	✓	✓	✓	✓	✓	✗	✓	✓	✓	✗	✗
Compression Range	0-75 mmHG	0-75 mmHG	0-50 mmHG	0-75 mmHg	0-50 mmHg	N/A	N/A	N/A	N/A	N/A	N/A
Sound Level Visual (Average highest level)					N/A	N/A	N/A	N/A	N/A	N/A	N/A
Sound Level dBa (Average highest level)	56 dBa	65 dBa	57 dBa	75 dBa	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Heat Therapy	✓	✗	✓	✗	✗	✗	✗	✗	✗	✗	✗
Heat Temperature	105°	N/A	105°	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Contrast Therapy	✗	✗	✓	✗	✗	✗	✗	✗	✗	✗	✗
Built in Timer	✓	✓	✓	✓	✗	✗	✗	✗	✗	✗	✓
Touch Screen	✓	✗	✓	✓	✗	✗	✗	✗	✗	✗	✓
DVT Sleeve Compatible	✓	✗	✓	✗	✗	✗	✗	✗	✗	✗	✗
Bilateral Wrap Compatible	✗	✓	✗	✗	✗	✓	✗	✗	✗	✗	✗